

pH Buffer Solutions



Choice of pH 4.01, pH 7.00 and pH 10.00 or custom formulations

Pint, Quart and Gallon Sizes

NIST Traceability Certificates Available at Additional Cost

Private/OEM Labeling Capability

Color-Coded for easy recognition

Sensorex supplies pH buffers as pre-mixed, ready-to-use, color-coded liquids.-- red (pH4.01), yellow (pH7.00) and blue (pH 10.00). All buffers are guaranteed accurate to within +/- 0.01 pH units @ 25 deg C. Available sizes include pints(473mL), quarts(946mL) and gallons(3.786L). Buffer labels include pH vs. temperature information for precise calibration. Custom buffers are also available on an OEM basis. Please contact Sensorex regarding your requirements.

Electrode storage solution is available in pint and gallon sizes. This is the same solution that is shipped with all pH electrodes and ensures electrodes are ready for immediate use. Storage solution is recommended for long-term electrode storage since it is not color-coded and will not impart any color to the electrode's reference junction.

Ordering Information

Description	Size
pH 4.01 +/- 0.01pH liquid buffer(RED)	Pint (473mL)
pH 4.01 +/- 0.01pH liquid buffer(RED)	Quart (946mL)
pH 4.01 +/- 0.01pH liquid buffer(RED)	US Gallon (3.785L)
pH 7.00 +/- 0.01pH liquid buffer(YELLOW)	Pint (473mL)
pH 7.00 +/- 0.01pH liquid buffer(YELLOW)	Quart (946mL)
pH 7.00 +/- 0.01pH liquid buffer(YELLOW)	US Gallon (3.785L)
pH 10.00 +/- 0.01pH liquid buffer(BLUE)	Pint (473mL)
pH 10.00 +/- 0.01pH liquid buffer(BLUE)	Quart (946mL)
pH 10.00 +/- 0.01pH liquid buffer(BLUE)	US Gallon (3.785L)
pH electrode storage solution(COLORLESS)	Pint (473mL)
pH electrode storage solution (COLORLESS)	US Gallon (3.785L)
	pH 4.01 +/- 0.01pH liquid buffer(RED) pH 4.01 +/- 0.01pH liquid buffer(RED) pH 4.01 +/- 0.01pH liquid buffer(RED) pH 7.00 +/- 0.01pH liquid buffer(YELLOW) pH 7.00 +/- 0.01pH liquid buffer(YELLOW) pH 7.00 +/- 0.01pH liquid buffer(YELLOW) pH 10.00 +/- 0.01pH liquid buffer(BLUE) pH 10.00 +/- 0.01pH liquid buffer(BLUE) pH 10.00 +/- 0.01pH liquid buffer(BLUE) pH electrode storage solution(COLORLESS)



Garden Grove, CA 92841 USA Tel: 714-895-4344 Fax: 714-894-4839 E-mail: info@sensorex.com www.sensorex.com